

AD-A118 334 ARMY MOBILITY EQUIPMENT RESEARCH AND DEVELOPMENT COMM--ETC F/G 7/4  
IDENTIFICATION OF DRYING OILS IN PAINTS BY GAS-LIQUID CHROMATOGRAPH--ETC(U)  
JUN 82 L. HERRIOTT  
UNCLASSIFIED MERADCOM-2360 NL

1 of 1  
AD-A  
118 334

END  
9-82

AD A118334

DTIC FILE COPY



AD

12


Report 2360

IDENTIFICATION OF DRYING OILS IN PAINTS  
BY GAS-LIQUID CHROMATOGRAPHY

by  
Laurie Herriott

June 1982

Approved for public release; distribution unlimited



U.S. ARMY MOBILITY EQUIPMENT  
RESEARCH AND DEVELOPMENT COMMAND  
FORT BELVOIR, VIRGINIA

82 08 17 063

**Destroy this report when it is no longer needed.  
Do not return it to the originator.**

---

**The citation in this report of trade names of  
commercially available products does not constitute  
official endorsement or approval of the use of such  
products.**

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

| REPORT DOCUMENTATION PAGE  |                                     | READ INSTRUCTIONS<br>BEFORE COMPLETING FORM                    |
|--|-------------------------------------|--|
| 1. REPORT NUMBER<br>2360   | 2. GOVT ACCESSION NO.<br>AD-A118334 | 3. RECIPIENT'S CATALOG NUMBER                                  |
| 4. TITLE (and Subtitle)<br>IDENTIFICATION OF DRYING OILS IN PAINTS BY<br>GAS-LIQUID CHROMATOGRAPHY   |                                     | 5. TYPE OF REPORT & PERIOD COVERED<br>Interim                  |
|  |                                     | 6. PERFORMING ORG. REPORT NUMBER                               |
| 7. AUTHOR(s)<br>Laurie Herriott  |                                     | 8. CONTRACT OR GRANT NUMBER(s)                                 |
| 9. PERFORMING ORGANIZATION NAME AND ADDRESS<br>Material Technology Lab; ATTN: DRDME-VO;<br>US Army Mobility Equipment Research and Development<br>Command; Fort Belvoir, VA 22060  |                                     | 10. PROGRAM ELEMENT, PROJECT, TASK<br>AREA & WORK UNIT NUMBERS |
| 11. CONTROLLING OFFICE NAME AND ADDRESS  |                                     | 12. REPORT DATE<br>June 1982                                   |
|  |                                     | 13. NUMBER OF PAGES<br>22                                      |
| 14. MONITORING AGENCY NAME & ADDRESS (If different from Controlling Office)  |                                     | 15. SECURITY CLASS. (of this report)<br>Unclassified           |
|  |                                     | 15a. DECLASSIFICATION/DOWNGRADING<br>SCHEDULE                  |
| 16. DISTRIBUTION STATEMENT (of this Report)<br><br>Approved for public release; distribution unlimited.  |                                     |  |
| 17. DISTRIBUTION STATEMENT (of the abstract entered in Block 20, if different from Report)   |                                     |  |
| 18. SUPPLEMENTARY NOTES  |                                     |  |
| 19. KEY WORDS (Continue on reverse side if necessary and identify by block number)   |                                     |  |
| 20. ABSTRACT (Continue on reverse side if necessary and identify by block number)<br><br>This gas-liquid chromatographic method is for the identification of drying oils in paint.<br>The identification is made from the distribution of fatty acids as methyl esters shown<br>on the chromatogram. |                                     |  |

DD FORM 1 JAN 73 1473


EDITION OF 1 NOV 65 IS OBSOLETE

UNCLASSIFIED

SECURITY CLASSIFICATION OF THIS PAGE (When Data Entered)

**UNCLASSIFIED**

**SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)**



**UNCLASSIFIED**

**ii**

**SECURITY CLASSIFICATION OF THIS PAGE(When Data Entered)**

## PREFACE

The work covered by this report was performed by Laurie Herriott and was reviewed by Troy Nichols under the supervision of Emil J. York, Chief, Material Technology Laboratory. The effort constitutes a part of the U.S. Army Materials and Mechanics Research Center's program to provide better procedures in chemical testing techniques to increase inspection efficiency for products procured by the U.S. Army Materiel Development and Readiness Command.



|                    |  |
|--------------------|--|
| Accession For      |  |
| NTIS GRA&I         | <input checked="checked" type="checkbox"/> |
| DTIC TAB           | <input type="checkbox"/>                   |
| Unannounced        | <input type="checkbox"/>                   |
| Justification      |  |
| By                 |  |
| Distribution/      |  |
| Availability Codes |  |
| Dist               | Avail and/or<br>Special                    |
| A                  |  |

## CONTENTS

| Section | Title                                | Page |
|---------|--------------------------------------|------|
|         | PREFACE                              | iii  |
|         | ILLUSTRATIONS                        | v    |
|         | TABLES                               | v    |
| I       | INTRODUCTION                         |      |
|         | 1. Purpose                           | 1    |
|         | 2. Background                        | 1    |
| II      | DETAILS OF TESTS                     |      |
|         | 3. Reagents                          | 1    |
|         | 4. Gas-Liquid Chromatographic Column | 1    |
|         | 5. Sample Preparation                | 2    |
|         | 6. Identification of Ester and Oil   | 3    |
| III     | RESULTS AND DISCUSSION               | 5    |
| IV      | CONCLUSION                           | 5    |

## ILLUSTRATIONS

| Figure | Title                                    | Page |
|--------|--|------|
| 1      | Chromatogram of Soybean Oil Methyl Ester | 6    |
| 2      | Chromatogram of Castor Oil Methyl Ester  | 7    |

## TABLES

| Table | Title   | Page |
|-------|---|------|
| 1     | Drying Oils Used in Sample Preparation                        | 2    |
| 2     | Relative Retention Data of Methyl Esters of Oil Fatty Acids   | 3    |
| 3     | Weight Percent Distribution of Fatty Acids in Unmodified Oils | 4    |



# IDENTIFICATION OF DRYING OILS IN PAINTS

## BY GAS-LIQUID CHROMATOGRAPHY

### I. INTRODUCTION

**1. Purpose.** The purpose of this investigation was to revise and update Federal Test Method Standard (FTMS) 7501, "Identification of Drying Oils (Gas-Liquid Chromatography)."

**2. Background.** Some of the military specifications for paints specify that certain drying oils are to be used in the formulation. If those drying oils are not used, the paint will prove unsatisfactory in its intended application. Therefore, it is essential for purposes of quality control that an accurate method be available for the identification of these drying oils. In September 1965, Method 7501 of FTMS 141 was approved for this purpose, but it has since been shown to be lacking in necessary flexibility for application to actual samples of specification paints.

### II. DETAILS OF TESTS

#### 3. Reagents.

a. Boron Trifluoride solution<sup>1</sup>—add 23 ml of boron fluoride methyl ether slowly to 150 ml of chilled absolute methyl alcohol in a 200-ml volumetric flask and dilute the volume with absolute methanol. Stopper tightly and store in a refrigerator. Discard after one month.

b. Petroleum Ether.<sup>2</sup>

**4. Gas-Liquid Chromatographic Column.** Prepare a 6-ft by 1/8-in. glass column, packed with 20-percent diethylene glycol succinate on 60-80 mesh chromosorb W acid washed (WAW). Using the thermal conductivity detector (TCD), condition the column at 200° C with helium at a flow rate of 60 ml/min, until the baseline is stabilized.

<sup>1</sup> L. D. Metcalfe and A. A. Schmitz, J. Am. Oil Chem. Soc., 33, 361 (1961).

<sup>2</sup> L. D. Metcalfe, A. A. Schmitz, and J. R. Felka, Anal. Chem., 38, 514 (1966).

**5. Sample Preparation.** Isolate the fatty acids from the resin or oil by saponification, acid hydrolysis, and solvent separation as in FTMS 7021 and FTMS 7031. If the fatty acids have become solid, warm them slightly. Then, accurately weigh about 0.2 g in a 6-in. screw cap test tube. Add 6 ml of the Boron Trifluoride solution and immerse in a 65° to 70° C water bath for 10 min. Transfer to a separatory funnel with 20 ml water and 30 ml of petroleum ether. Shake gently. Discard the bottom water layer and filter the remaining layer through fast filter paper into a 50-ml beaker and evaporate in a 65° to 70° C water bath to about 2 ml. Do not evaporate to dryness. Transfer and store in a closed 2-ml vial. Monitoring the specification or statement of composition, inject 1.5 ml of the prepared methyl esters using the optimum conditions for each oil as shown in Table 1.

Table 1. Drying Oils Used in Sample Preparation

| Drying Oils | Column Temperature<br>(° C)                          | Carrier (Helium) Flow<br>(ml/min) |
|-------------|--|-----------------------------------|
| Castor      | 20 min at 170°;<br>50 min at 200°;<br>rate 20° C/min | 60                                |
| Soya        | 180  | 60                                |
| Cottonseed  | 180  | 60                                |
| Linseed     | 180  | 50                                |
| Safflower   | 190  | 40                                |

Thermal conductivity detector 250° C.  
Injection port temperature 275° C.

**6. Identification of Ester and Oil.** The individual methyl esters can be identified by their relative retention times as given in Table 2. The drying oil, itself, may be identified in either of two ways. First, by comparison of its chromatogram with those of known oils. Second, identification of the drying oil may be made by matching concentration of the methyl esters with the values given in Table 3. For this purpose, concentration ( $C_x$ ) of the methyl ester is calculated as follows:

$$C_x = \frac{A_x}{A_T} \times F_x \times 100.$$

Where  $A_x$  is the peak area of the methyl ester;  $A_T$  is the total area of the methyl ester peaks; and  $F_x$  is a correction factor determined from a known sample.

Table 2. Relative Retention Data of Methyl Esters of Oil Fatty Acids

| Ester       | Relative Retention Time |
|-------------|-------------------------|
| Caprylic    | 0.10                    |
| Capric      | 0.19                    |
| Lauric      | 0.32                    |
| Myristic    | 0.56                    |
| Palmitic    | 1.0                     |
| Palmitoleic | 1.15                    |
| Stearic     | 1.7                     |
| Oleic       | 1.94                    |
| Linoleic    | 2.31                    |
| Linolenic   | 2.97                    |
| Licanic     | 5.2                     |
| Eleostearic | 5.3                     |
| Ricinoleic  | 10.0                    |

(Ester of Palmitic Acid = 1)

Table 3. Weight Percent Distribution of Fatty Acids in Unmodified Oil

| Oil        | Caprylic | Capric | Lauric | Myristic | Palmitic | Palmitoleic | Stearic | Oleic | Linoleic | Linolenic | Licanic | Eleostearic | Ricinoleic |
|------------|----------|--------|--------|----------|----------|-------------|---------|-------|----------|-----------|---------|-------------|------------|
| Castor     |          |        |        |          | 2        |             | 1       | 7     | 3        |           |         |             | 87         |
| Coconut    | 6        | 6      | 44     | 18       | 11       |             | 6       | 7     | 2        | t         |         |             |            |
| Cottonseed |          | t      | t      | 1        | 29       | 2           | 4       | 24    | 40       |           |         |             |            |
| Linseed    |          |        |        |          | 6        | t           | 4       | 22    | 16       | 52        |         |             |            |
| Oiticica   |          |        |        |          | 7        |             | 5       | 6     |          |           |         |             | 78         |
| Perilla    |          |        |        |          | 7        |             | 2       | 13    | 14       | 64        |         |             |            |
| Safflower  |          |        |        | t        | 8        | t           | 3       | 13    | 75       | 1         |         |             |            |
| Soybean    |          |        |        | t        | 11       |             | 4       | 25    | 51       | 9         |         |             |            |
| Tall       |          |        |        |          | 5        |             | 3       | 46    | 41       | 3         |         |             |            |
| Tung       |          |        |        |          | 4        |             | 1       | 8     | 4        | 3         |         |             |            |
| Menhaden   |          | t      | 7      | 16       | 16       | 2           | 15      | 7     | 2        |           |         |             |            |

t = trace

### III. RESULTS AND DISCUSSION

The various drying oils used in paints are composed of specific amounts of fatty acids as given in Table 3. Determination of these fatty acids as methyl esters by gas-liquid chromatography provides a means for confirmation of the required drying oil in the paint.

Because of the wide variation of composition of drying oils it was impractical, if not impossible, to develop a single chromatographic program for satisfactory identification for all of the oil acids. Carrier flow rate and column oven temperature were adjusted for optimum resolution of the methyl esters for each drying oil.

The one complete method for identification of the oil acid consists of the following steps. The acid is saponified and extracted with water as the sodium salt from an ethyl ether despersion of the paint vehicle. On acidification with hydrochloric acid, the fatty acid is then extracted with chloroform (FTMS 7021). Evaporation<sup>1 2</sup> follows with conversion of the fatty acid to the methyl ester with Boron Trifluoride solution. This well-known reaction converts the acid to a more volatile product which can then be determined by gas-liquid chromatography. An example of this reaction is:




Figure 1 is a chromatogram of the methyl esters of soybean oil obtained with the above procedure. This chromatogram illustrates the use of this method in the identification of drying oils. By its retention time the methyl ester is identified, and by the percentage of each methyl ester, the drying oil is identified.


Figure 2 is the chromatogram of the methyl esters of castor oil. The striking difference between the chromatograms of Figures 1 and 2 illustrate the utilization of this method in identifying drying oils.

### V. CONCLUSION

An improved method was developed for the confirmation of the presence of the required drying oil in military specification paints. To a limited extent, the method is useful in the identification of drying oils in nonspecification paints. In some instances, mixtures of drying oils can be determined.


<sup>1</sup> L. D. Metcalfe and A. A. Schmitz, J. Am. Oil Chem. Soc., 33, 361 (1961).

<sup>2</sup> L. D. Metcalfe, A. A. Schmitz, and J. R. Pelka, Anal. Chem., 38, 514 (1966).



1. Petroleum Ether
2. Palmitic
3. Stearic
4. Oleic
5. Linoleic
6. Linolenic

Figure 1. Chromatogram of soybean oil methyl ester.



1. Petroleum Ether
2. Palmitic
3. Stearic
4. Oleic
5. Linoleic
6. Impurity
7. Ricinoleic

Figure 2. Chromatogram of castor oil methyl ester.

# **DISTRIBUTION FOR MERADCOM REPORT 2360**

| <b>No. Copies</b> | <b>Addressee</b>   | <b>No. Copies</b> | <b>Addressee</b>  |
|-------------------|--|-------------------|---|
| 2                 | Metals and Ceramics Information Center<br>ATTN: Mr. Harold Mindlin<br>Director<br>Mr. James Lynch<br>Assistant Director<br>505 King Avenue<br>Columbus, OH 43201 | 3                 | Commander<br>US Army Electronics R&D Command<br>ATTN: DRSEL-PA-E<br>Mr. Stan Alster<br>Mr. J. Quinn<br>DELS-D<br>Fort Monmouth, NJ 07703  |
| 12                | Defense Technical Information Center<br>Cameron Station<br>Alexandria, VA 22314  | 9                 | Commander<br>US Army Missile Research and Development Command<br>ATTN: DRDMI-TB, Redstone Scientific Info Ctr (2)<br>DRDMI-TK, Mr. J. Alley<br>DRDMI-M<br>DRDMI-ET, Mr. Robert O. Black<br>DRDMI-QS, Mr. George L. Stewart, Jr.<br>DRDMI-EAT, Mr. R. Talley<br>DRDMI-QP<br>DRSMI-RR<br>Redstone Arsenal, AL 35809 |
| 1                 | Commander<br>US Army Foreign Science and Technology Center<br>ATTN: DRXST-SD3<br>220 Seventh Street, NE<br>Charlottesville, VA 22901                             | 4                 | Commander<br>US Army Troop Support and Aviation Materiel Readiness Command<br>ATTN: DRSTS-PL (2)<br>Mr. J. Corwin (1)<br>DRSTS-Q<br>DRSTS-MES (1)<br>DRSTS-M<br>4300 Goodfellow Blvd<br>St. Louis, MO 63120   |
| 2                 | Office of the Deputy Chief of Staff for Research, Development, and Acquisition<br>ATTN: DAMA-ARZ-E<br>DAMA-CSS<br>Washington, DC 20310                           | 1                 | Commander<br>US Army Natick Research and Development Command<br>ATTN: DRDNA-EM<br>Natick, MA 01760  |
| 2                 | Commander<br>Army Research Office<br>ATTN: Dr. George Mayer<br>Mr. J. J. Murray<br>P. O. Box 12211<br>Research Triangle Park, NC 27709                           | 2                 | Commander<br>US Army Tank-Automotive Materiel Readiness Command<br>ATTN: DRSTA-Q<br>Warren, MI 48090  |
| 6                 | Commander<br>US Army Materiel Development and Readiness Command<br>ATTN: DRCQA-E<br>DRCDE-D<br>DRCDMD-FT<br>DRCLDC<br>DRDCMT<br>DRCMM-M<br>Alexandria, VA 22333  |                   |   |



| No. Copies | Addressee  | No. Copies | Addressee   |
|------------|--|------------|---|
| 6          | Commander<br>US Army Armament Materiel<br>Readiness Command<br>ATTN: DRSAR-QA (2)<br>DRSAR-SC<br>DRSAR-RDP<br>DRSAR-EN<br>DRSAE-QAE<br>Rock Island, IL 61299   | 1          | Commander<br>US Army Aberdeen Proving Ground<br>ATTN: STEAP-MT-U (GE Branch)<br>Aberdeen Proving Ground, MD<br>21005        |
| 4          | Commander<br>Rock Island Arsenal<br>ATTN: SARRI-EN, Mr. W. M.<br>Kisner<br>SARRI-ENM, Mr. W. D.<br>McHenry<br>SARRI-QA<br>SARRI-LPL<br>Rock Island, IL 61299   | 1          | Director<br>US Army Materiel Systems Analysis<br>Agency<br>ATTN: DRXSY-CM<br>Aberdeen Proving Ground, MD<br>21005           |
| 12         | Commander<br>US Army Armament Research<br>and Development Command<br>ATTN: DRDAR-LC, Mr. E. Kelly<br>DRDAR-LCA, Dr. Sharkoff<br>DRDAR-LCE, Dr. Walker<br>DRDAR-QAS, Mr. F.<br>Fitzsimmons (5)<br>DRDAR-SCM, Jr. J. Corrie<br>DRDAR-TSP, Mr. B.<br>Stephans<br>DRDAR-TSS (STINFO)(2)<br>Dover, NJ 07703 | 1          | Director<br>US Army Ballistic Research Lab<br>ATTN: DRDAR-TSD-S (STINFO)<br>Aberdeen Proving Ground, MD 21005               |
| 2          | Commander<br>Edgewood Arsenal<br>ATTN: DRDAR-CLR, Mr.<br>Montanary<br>DRDAR-QAC, Dr. Maurits<br>Aberdeen Proving Ground, MD<br>21010   | 2          | Commander<br>Watervliet Arsenal<br>ATTN: DRDAR-LCB<br>Mr. T. Moraczewski<br>SARWV-PPI, Mr. L. Jette<br>Watervliet, NY 12189 |
| 1          | Technical Library<br>Chemical Systems Laboratory<br>Aberdeen Proving Ground, MD<br>21010   | 4          | Commander<br>US Army Aviation R&D Command<br>ATTN: DRDAV-EXT<br>DRDAV-QR<br>DRDAV-QP<br>DRDAV-QE<br>St. Louis, MO 63166     |

| No. Copies | Addressee   | No. Copies | Addressee   |
|------------|---|------------|---|
| 7          | Commander<br>US Army Tank-Automotive<br>Research and Development Cnd<br>ATTN: DRDTA-UL, Tech Lib<br>DRDTA-RCKM, Mr. S.<br>Goodman<br>DRDTA-RCKT, Mr. J. Fix<br>DRDTA-RTAS, Mr. S.<br>Catalano<br>DRDTA-TTM, Mr. W.<br>Moncrief<br>DRDTA-ZA, Mr. O. Renius<br>(Only Infrared, Ultrasonic,<br>or Holographic Reports)<br>DRDTA-JA, Mr. C. Kedzior<br>Warren, MI 48090 | 1          | Commander<br>US Army Yuma Proving Ground<br>ATTN: Technical Library<br>Yuma, AZ 85364   |
|            |   | 1          | Commander<br>US Army Tropic Test Center<br>ATTN: STETC-TD, Drawer 942<br>Fort Clayton, Canal Zone   |
|            |   | 4          | Commander<br>Aberdeen Proving Ground<br>ATTN: STEAP-MT<br>STEAP-TL<br>STEAP-MT-M, Mr. J. A.<br>Feroi<br>STEAP-MT-G, Mr. R. L.<br>Huddleston<br>Aberdeen Proving Ground, MD<br>21005 |
| 1          | Director<br>US Army Industrial Base Engineering<br>Activity<br>ATTN: DRXIB-MT, Mr. D. Brim<br>Rock Island, IL 61299   | 1          | Commander<br>US Army Cold Region Test Center<br>ATTN: STECR-OP-PM<br>APO Seattle, Washington 98733  |
| 1          | Commander<br>Harry Diamond Laboratories<br>ATTN: DELHD-EDE<br>Mr. B. F. Willis<br>2800 Powder Mill Road<br>Adelphi, MD 20783  | 1          | Commander<br>US Army Dugway Proving Ground<br>ATTN: STEDP-MT<br>Dugway, UT 84022  |
| 2          | Commander<br>US Army Test and Evaluation<br>Command<br>ATTN: DRSTE-TD<br>DRSTE-ME<br>Aberdeen Proving Ground, MD<br>21005   | 1          | Commander<br>US Army Electronic Proving Ground<br>ATTN: STEEP-MT<br>Fort Huachuca, AZ 85613   |
| 4          | Commander<br>US Army White Sands Missile<br>Range<br>ATTN: STEWS-AD-L<br>STEWS-ID<br>STEWS-TD-PM<br>ATAA-SL (Tech Lib)<br>White Sands Missile Range, NM<br>88002  | 1          | Commander<br>Jefferson Proving Ground<br>ATTN: STEJP-TD-I<br>Madison, IN 47250  |
|            |   | 1          | Commander<br>US Army Aircraft Development<br>Test Activity<br>ATTN: STEBG-TD<br>Fort Rucker, AL 36362   |

| No. Copies | Addressee   | No. Copies | Addressee  |
|------------|---|------------|--|
| 2          | President<br>US Army Armor and Engineer<br>Board<br>ATTN: ATZKOE-TA<br>ATZK-AE-PD-E<br>Fort Knox, KY 40121            | 1          | Commander<br>Sacramento Army Depot<br>ATTN: SDSSA-QA<br>Sacramento, CA 95813                                 |
| 1          | President<br>US Army Field Artillery Board<br>ATTN: ATZA-BDOP<br>Fort Sill, OK 73503                                  | 1          | Commander<br>Savanna Army Depot Activity<br>ATTN: SDSSV-S<br>Savanna, IL 61074                               |
| 1          | Commander<br>Anniston Army Depot<br>ATTN: SDSAN-QA<br>Anniston, AL 36202  | 1          | Commander<br>Seneca Army Depot<br>ATTN: SDSEE-R<br>Romulus, NY 14541   |
| 1          | Commander<br>Corpus Christi Army Depot<br>ATTN: SDSCC-MEE<br>Mr. Haggerty<br>Mail Stop 55<br>Corpus Christi, TX 78419 | 1          | Commander<br>Sharpe Army Depot<br>ATTN: SDSSH-QE<br>Lathrop, CA 95330  |
| 1          | Commander<br>Letterkenny Army Depot<br>ATTN: SDSLE-QA<br>Chambersburg, PA 17201                                       | 1          | Commander<br>Sierra Army Depot<br>ATTN: SDSSI-DQA<br>Herlong, CA 96113                                       |
| 1          | Commander<br>Lexington-Bluegrass Army Depot<br>ATTN: SDSRR-QA<br>Lexington, KY 40507                                  | 1          | Commander<br>Tobyhanna Army Depot<br>ATTN: SDSTO-Q<br>Tobyhanna, PA 18466                                    |
| 1          | Commander<br>New Cumberland Army Depot<br>ATTN: SDSNC-QA<br>New Cumberland, PA 17070                                  | 1          | Commander<br>Tooele Army Depot<br>ATTN: SDSTE-QA<br>Tooele, UT 84074   |
| 2          | Commander<br>US Army Depot Activity<br>Pueblo<br>ATTN: SDSTE-PU-Q<br>Pueblo, CO 81001                                 | 1          | Director<br>DARCOM Ammunition Center<br>ATTN: SARAC-DC<br>Savanna, IL 61704                                  |
| 1          | Commander<br>Red River Army Depot<br>ATTN: SDSRR-QA<br>Texarkana, TX 75501  | 2          | Naval Research Laboratory<br>ATTN: Dr. J. M. Krafft, Code 8430<br>Library, Code 2620<br>Washington, DC 20375 |

| No. Copies | Addressee  | No. Copies | Addressee  |
|------------|--|------------|--|
| 3          | Air Force Materials Laboratory<br>ATTN: AFML-DO, Library<br>AFML-LTM, Mr. E. Wheeler<br>AFML-LLP, Mr. R. Rowand<br>Wright-Patterson AFB, OH 45433  | 1          | HQDA (DAEN-RDL)<br>Washington, DC 20314  |
|            |  | 1          | HQDA (DAEN-MPE-T)<br>Washington, DC 20314  |
| 9          | Director<br>Army Materials and Mechanics<br>Research Center<br>ATTN: DRXMR-P (1)<br>DRXMR-PL (2)<br>DRXMR-M (2)<br>DRXMR-MQ (1)<br>DRXMR-MI, Mr. Darcy (1)<br>DRDMM-L, Dr. Chait (1)<br>DRXMR-RA, Mr. Valente (1)<br>DRXMR-AG-MD (1)<br>DRXMR-X (1)<br>DRXMR-PR (1)<br>DRXMR-T (1)<br>DRXMR-E (1)<br>Watertown, MA 02172 | 1          | Director<br>US Army Engineer Waterways Experiment<br>Station<br>ATTN: Chief, Library Branch<br>Technical Information Center<br>Vicksburg, MS 39180 |
|            |  | 2          | Director<br>Petrol & Fld Svc Dept<br>US Army Quartermaster School<br>Fort Lee, VA 23801  |
|            |  | 1          | President<br>US Army Aviation Test Board<br>ATTN: STEBG-PO<br>Fort Rucker, AL 36360  |
| 1          | Director, Technical Information<br>Defense Advanced Research Projects<br>Agency<br>1400 Wilson Blvd<br>Arlington, VA 22209   | 1          | US Army Aviation School Library<br>P. O. Drawer O<br>Fort Rucker, AL 36360   |
|            |  | 2          | HQ, 193D Infantry Brigade (Pan)<br>ATTN: AFZU-FE<br>APO Miami 34004  |
| 1          | Director<br>Defense Nuclear Agency<br>ATTN: TITL<br>Washington, DC 20305   | 2          | Special Forces Detachment, Europe<br>ATTN: PBO<br>APO New York 09050   |
| 1          | Commander, HQ TRADOC<br>ATTN: ATEN-ME<br>Fort Monroe, VA 23651   | 2          | Engineer Representative<br>USA Research & Standardization<br>Group (Europe)<br>Box 65<br>FPO 09510   |
| 1          | HQDA (DAMA-AOA-M)<br>Washington, DC 20310  |            |  |
| 1          | HQDA (DALO-TSM)<br>Washington, DC 20310  | 1          | HQDA<br>ODCSLOG<br>DALO-TSE<br>Room 1E588<br>Pentagon, Washington, DC 20310  |

| No. Copies | Addressee   | No. Copies | Addressee   |
|------------|---|------------|---|
| 1          | Plastics Technical Evaluation Center<br>ARRADCOM, Bldg 3401<br>ATTN: A. M. Anzalone<br>Dover, NJ 07801              |            | Spec Asst/Matl Asmt, DRDME-ZG<br>Spec Asst/Scs & Tech, DRDME-ZK<br>CIRCULATE  |
| 1          | Commander<br>Frankford Arsenal<br>ATTN: Library, K2400, B151-2<br>Philadelphia, PA 19137                            | 1          | C, Ctrmine Lab, DRDME-N<br>C, Engy & Wtr Res Lab, DRDME-G<br>C, Elec Pwr Lab, DRDME-E<br>C, Camo & Topo Lab, DRDME-R<br>C, CS/CI Lab, DRDME-M<br>C, Mech & Constr Eqpt Lab, DRDME-H<br>C, Ctr Intrus Lab, DRDME-X<br>C, Matl Tech Lab, DRDME-V<br>Dir, Prod A&T Dir, DRDME-T<br>CIRCULATE |
| 1          | Commandant<br>US Army Engineer School<br>ATTN: ATZA-CDD<br>Fort Belvoir, VA 22060                                   |            |   |
| 1          | President<br>US Army Airborne, Communications<br>and Electronics<br>ATTN: STEBF-ABTD<br>Fort Bragg, NC 28307        | 5          | Material Technology Lab, DRDME-V  |
|            |   | 10         | Organic Chem Ctg Rsch Grp, DRDME-VO   |
|            |   | 3          | Tech Reports Ofc, DRDME-WP  |
|            |   | 3          | Security Ofc (for liaison officers),<br>DRDME-S   |
|            |   | 2          | Tech Library, DRDME-WC  |
| 1          | Commander<br>Headquarters, 39th Engineer<br>Battalion (Cbt)<br>Fort Devens, MA 01433                                | 1          | Programs & Anal Dir, DRDME-U  |
|            |   | 1          | Pub Affairs Ofc, DRDME-I  |
|            |   | 1          | Office of Chief Counsel, DRDME-L  |
| 1          | Commander and Director<br>USA FESA<br>ATTN: FESA-TS<br>Fort Belvoir, VA 22060                                       | 1          | Director, Physics Program (421)<br>Office of Naval Research<br>Arlington, VA 22217  |
|            |   | 2          | Commander, Naval Facilities<br>Engineering Command<br>Department of the Navy<br>ATTN: Code 032-B, 062<br>200 Stovall St<br>Alexandria, VA 22332   |
| 1          | HQ, USAEUR & Seventh Army<br>Deputy Chief of Staff, Engineer<br>ATTN: AEAEN-MT-P<br>APO New York 09403              |            |   |
| 1          | HQ, USAEUR & Seventh Army<br>Deputy Chief of Staff, Operations<br>ATTN: AEAGC-FMD<br>APO New York 09403             | 1          | US Naval Oceanographic Office<br>Navy Library/NSTL Station<br>Bay St. Louis, MS 39522   |
|            | MERADCOM  | 1          | Library (Code L08A)<br>Civil Engineering Laboratory<br>Naval Construction Battalion Ctr<br>Port Hueneme, CA 93043   |
| 1          | Commander, DRDME-Z<br>Tech Dir, DRDME-ZT<br>Assoc Tech Dir/R&D, DRDME-ZN<br>Assoc Tech Dir/Engrg & Acq,<br>DRDME-ZE |            |   |

| No. Copies | Addressee  |
|------------|--|
| 1          | Director<br>Earth Physics Program<br>Code 464<br>Office of Naval Research<br>Arlington, VA 22217                                 |
| 1          | Naval Training Equipment Center<br>ATTN: Tech Library<br>Orlando, FL 32813   |
| 1          | HQ USAF/RDPT<br>ATTN: Mr. Allan Eaffy<br>Washington, DC 20330  |
| 1          | HQ USAF/LEEEU<br>Chief, Utilities Branch<br>Washington, DC 20330   |
| 1          | US Air Force<br>HQ Air Force Engineering and<br>Services Center<br>Technical Library FL 7050<br>Tyndall AFB, FL 32403            |
| 1          | Chief, Lubrication Br<br>Fuels & Lubrication Div<br>ATTN: AFWAL/POSL<br>Wright-Patterson AFB, OH 45433                           |
| 1          | Department of Transportation<br>Library, FOP 10A, M494-6<br>800 Independence Ave, SW<br>Washington, DC 20591                     |
| 1          | Professor Raymond R. Fox<br>School of Engineering and Applied<br>Science<br>George Washington University<br>Washington, DC 20052 |
| 1          | Reliability Analysis Center<br>Rome Air Development Center<br>RADC/RBRAC (I. L. Krulac)<br>Griffiss AFB, NY 13441                |